

Genomic DNA Isolation from Whole Blood Kit (KGI)

Introduction

Alpha Biobeads' blood DNA kit is designed to isolate nucleic acid from human whole blood. The magnetic nanoparticle (1 μm) technology enables efficient, consistent, and scalable extraction of genomic DNA from as little as 10 μL and up to 10mL blood. The obtained DNA is with high yield and purity (isolate 10 to 20 μg DNA from 350 μl blood, A260/280 >1.7, A260/230 >1.8), and is suitable for direct downstream applications such as restriction enzyme digestion, PCR, qPCR and NGS. The magnetic particles feature uniformed size and slow sedimentation rate, which makes it an excellent choice for automation system applicable to high throughput sample preparation.

Kit Components and Storage


Each kit contains reagents for 96 reactions (based on 350 μl blood/reaction):

Kit Components	Quantity	Storage
Magnetic Beads	2.5 mL (40 mg beads/mL)	2 to 8 °C, do not freeze
Lysis Buffer (GLB-40)	40 mL	2 to 8 °C
Washing Buffer 1 (GWB-1)	2 x 60 mL	2 to 8 °C
Washing Buffer 2 (GWB-2)	2 x 10 mL	2 to 8 °C
Elution Buffer (GEB-15)	15 mL	2 to 8 °C

Materials to be Supplied by the User

- Human Whole Blood (Fresh or Frozen)
- Proteinase K
- Isopropanol or 2-Propanol
- Ethanol
- Magnetic Separator

Protocol Flow Chart



Critical Notes Before You Start:

- Reconstitute Washing Buffer 1: Add 40 mL Isopropanol to the Washing Buffer 1 before your first use. Mix well and store at 2 to 8 °C.
- Reconstitute Washing Buffer 2: Add 90 mL Ethanol to the Washing Buffer 2 before your first use. Mix well and store at 2 to 8 °C.
- Prepare Proteinase K Solution (20 mg/mL): Dissolve Proteinase K powder in 10 mM Tris-HCl buffer (pH 8.0).
- The current protocol is for genomic DNA isolation from 350 μL human whole blood. The protocol has been validated to scale up or down to different blood volumes. If the amount of blood used has been changed, the amount of magnetic beads, proteinase K solution, lysis buffer, washing buffers, and elution buffer will need to be adjusted accordingly.
- There might be beads retained in the pipette tip during pipetting steps. Carefully pipette a few more times to release the beads back into the reaction in order to maximize yield.
- For any vortex steps, vortex at maximum speed to ensure mixing.

Protocol

1. Add 50 μ L Proteinase K solution (20 mg/mL) to a 350 μ L blood sample in a 1.5 mL microcentrifuge tube. Vortex for 30 seconds and incubate for 2 minutes at room temperature.
2. Add 350 μ L lysis buffer to the sample, and then vortex for 30 seconds to mix.
3. Incubate at 55°C for 5 minutes, and then vortex for 30 seconds to mix.
4. Repeat step #3 one more time.
5. Add 25 μ L magnetic beads to the sample, and vortex for 15 seconds.
Note: Please re-suspend the magnetic beads before adding the beads to the blood sample.
6. Add 500 μ L isopropanol to the mixture, and vortex for 30 seconds to mix.
7. Incubate the mixture on a roller for 5 minutes, and vortex for 30 seconds to mix.
8. Place the tube on the magnetic separator. After 2 minutes, remove the supernatant using a pipette.
Note: Please try to avoid touching the magnetic beads pellet with the pipette.
9. Remove the tube from the magnetic separator and add 1 mL Washing Buffer 1 to the pellet. Pipette up and down to fully resuspend the beads and then vortex for 30 seconds to mix.
10. Repeat steps #8 and #9 one more time.
11. Place the tube on the magnetic separator. After 2 minutes, remove the supernatant using a pipette.
12. Remove the tube from the magnetic separator and add 1 mL Washing Buffer 2 to the pellet. Pipette up and down to fully resuspend the beads. Vortex for 30 seconds to mix.
13. Use the pipette to transfer the magnetic beads to a clean tube.
14. Place the tube on a magnetic separator for 1 minute. Remove the supernatant. Avoid touching the pellet with the pipette.
15. Remove the tube from the magnetic separator and add 1 mL Washing Buffer 2 and vortex for 30 seconds. Place the tube on the magnetic separator for 1 minute. Remove the supernatant.
16. Leave the tube on the magnetic separator for 5 minutes to air dry.
Note: Please remove any visible supernatant left in the tube with a small pipette tip. Don't let the pellet dry for over 5 minutes.
17. Take out the tube and add 100 μ L Elution Buffer to the pellet and vortex for 2 minutes.
Note: If it's difficult to resuspend the magnetic beads, a pipette could be used to pipette up and down to fully resuspend the beads.
18. Place the tube on the magnetic separator. After one minute, transfer the supernatant to a clean tube. This supernatant contains the purified genomic DNA.

WARRANTIES AND DISCLAIMER

The Alpha Biobeads products ("Product") is warranted to operate or perform in conformance with published Product specifications at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products **are supplied for research use only**. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non-infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

© 2016 Alpha Biobeads, LLC. All rights reserved. Unless otherwise indicated, all trademarks are property of Alpha Biobeads, LLC. Printed in the USA.